

DEFA6-defensin, alpha 6, Paneth cell-specific Gene

Reference: PTXBC069667

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DEF6; HD-6; defensin-6; defensin, alpha 6, Paneth cell-specific; defensin alpha 6		

For research use only.

GENE DESCRIPTION	defensin, alpha 6, Paneth cell-specific
GENE ID	1671
ORIGIN SPECIES	Human
ACCESSIONS	BC069667
VECTOR	pDONR223
NCBI SYMBOL	DEFA6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagaacccctaccatcctcactgctgttctctcgtggccctccaggccaaggctgagccactccaagctgaggatgatccactgcaggcaaaagcttataggctgatgccaggagcagcgtggggcaaatgaccaggactttgccgtctcctttgcagaggatgcaagctcaagtc
ttagagctttgggtcaacaagggttccacttgcagaaggtcctgtattcaacagaatattctatgggacctgcactgtcatgggtattaaccacagattctgctgcctctga