

DEFA5-defensin, alpha 5, Paneth cell-specific Gene

Reference: PTXBC069690

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DEF5; HD-5; defensin-5; HD5(20-94); defensin, alpha 5, Paneth cell-specific; defensin, alpha 5, preproprotein; defensin alpha 5		

For research use only.

GENE DESCRIPTION	defensin, alpha 5, Paneth cell-specific
GENE ID	1670
ORIGIN SPECIES	Human
ACCESSIONS	BC069690
VECTOR	pDONR223
NCBI SYMBOL	DEFA5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggaccatgcacatctgtgcccattctctgtggccctgcaggcccgctgagtcactccaggaaagagctgatgaggtacaaccagaagcagctctggggaagacaaccaggacctgctatctcttgcaggaaatggactctctgctcttagaacctcaggtctcaggaagagccacctgctattgccgaactggccgtgtgtctacccgtgagtcctctccgggtgtgtgaatcagtgccgcctctacagactctgtgtcgtga