

CCL18-chemokine (C-C motif) ligand 18 (pulmonary and activation-regulated) Gene

Reference: PTXBC069700

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AMAC-1; AMAC1; CKb7; DC-CK1; DCCK1; MIP-4; PARC; SCYA18; C-C motif chemokine 18; CC chemokine PARC; CC chemokine ligand 18; alternative macrophage activation-associated CC chemokine 1; chemokine (C-C motif) ligand 18 (pulmonary and activation-regulated);		

For research use only.

GENE DESCRIPTION	chemokine (C-C motif) ligand 18 (pulmonary and activation-regulated)
GENE ID	6362
ORIGIN SPECIES	Human
ACCESSIONS	BC069700
VECTOR	pDONR223
NCBI SYMBOL	CCL18
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaggccttgacgtgccctcctgtcctcgtcaccatggccctcgtcctcgtgcacaagttggtaccaacaagagctcgtcgtcctcgtctatcctcctggcagattccaaaaagttcatagtgactattctgaaccagccccagtgccccaagccaggtgtcatcctccta
accaagagaggccggcagatctgtcgtgacccaataagaagtggtccagaaatacatcagcgacctgaagctgaatgcctga