

## CNR2-cannabinoid receptor 2 (macrophage) Gene

Reference: PTXBC069722

|                |                                                                                                          |                     |                |
|----------------|----------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                      | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                     | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | testis-dominant CNR2 isoform CB2; CB-2; CX5; cannabinoid receptor 2; cannabinoid receptor 2 (macrophage) |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | cannabinoid receptor 2 (macrophage)                                                               |
| <b>GENE ID</b>          | 1269                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC069722                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | CNR2                                                                                              |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggaggagctggtgacagagatgccaatggctccaaggatggctggattccaacctatgaaggattacatgacctgagtggtgtccccaagacagctgttctgtgtgtgacactcttctgggctgctaaagtgcctggagaacgtggctgtgctctatctgacctgtctccca  
ccaactccgcccgaagccctacactgttcatggcagcttggctgggctgacttctgccaagtgtgtcttttcatgacgcttggtaattccatgttttccatgggtgtggaattccaagctgtcttctgctgaagattggcagctgactatgaccttcacagcctctgtggtagcctct  
gtgacccgcaatgacagatacctctgctgctatccaccttctacaaagctctgtcaccctgtggaaggacactgtgacccctgggcatcatgtgggtctctcagcactagctcctacctgcccctcatgggatggactgtgtccaggccctgctctgagcttttccactgatcc  
ccaatgactacctgctgagctggctctgttcatgcccctctcttccggaatcatctacacatggtgcatgttctctggaaggcccatcagcatgtggccagctgtgtggtccaccaggacagggcaggtgccagggaatggcccgaatgaggtggatgtgaggtggccaagacctag  
ggctagtgttggctgtgctctcatctgtgttccagctgtgcccctcatggccacagcctggccactacgctcagtgaccaggtcaagaaggccttcttctgctccatgtgtgctcctcaactccatggtaaccctgtcatctatgctctacggagtgagagatccactcctctgc  
ccatcactgctgctcactggaagaagtgtgtgagggccttgggtcagaggcaaaagaagaagccccgagatcctcagtcaccgagacagagctgatgggaaatcactccgtggccagattccagagatctagacctctgtattgctga