

GAST-gastrin Gene

Reference: PTXBC069724

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GAS; preprogastrin		

For research use only.

GENE DESCRIPTION	gastrin
GENE ID	2520
ORIGIN SPECIES	Human
ACCESSIONS	BC069724
VECTOR	pDONR223
NCBI SYMBOL	GAST
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagcgaactatgtgtatgtgctgatctttgcactggctctggccgcttctctgaagcttcttgaagccccgctccagcagccagatgcacccttaggtacaggggccaacaggacgtggagctaccctggctggagcagcagggcccgctctcatcgaaggcagctgg
gaccccgagggtcccccacacctcgtggcagaccgtccaagaagcagggacccatggctggaggagaagaagaagcctatggatggatggactcggcccgccagtgctgaggatgagaactaa