

DSCR4-Down syndrome critical region gene 4 Gene

Reference: PTXBC069729

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DCRB; DSCRB; Down syndrome critical region protein 4; Down syndrome critical region gene 4; Down syndrome critical region protein B; Down syndrome critical region 4		

For research use only.

GENE DESCRIPTION	Down syndrome critical region gene 4
GENE ID	10281
ORIGIN SPECIES	Human
ACCESSIONS	BC069729
VECTOR	pDONR223
NCBI SYMBOL	DSCR4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgttaatcatcttgacgagagatgatgaaccccgatattaccacagacagtgatgccgttcaccagcattgcactctactccccctctgatcctgcctcagcttctctctccacagagaagaaaaattctgctaaagtctgcaacatggttctgctgagtttcagcctgcc
agcttctctacggattctgactfgccagccccacaataaccagagagggcagcaattttggccaaaatgtctgatttggaaataccaacttactccatgggctccacaagaatcagatgggagaagggacaagcagataagcgaagcccatcaacctga