

RPA4-replication protein A4, 34kDa Gene

Reference: PTXBC069791

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	replication protein A complex 34 kd subunit homolog Rpa4; HSU24186; replication protein A 30 kDa subunit; RF-A protein 4; RP-A p30; replication factor A protein 4; replication protein A4, 30kDa; replication protein A4, 34kDa; replication protein A4		

For research use only.

GENE DESCRIPTION	replication protein A4, 34kDa
GENE ID	29935
ORIGIN SPECIES	Human
ACCESSIONS	BC069791
VECTOR	pDONR223
NCBI SYMBOL	RPA4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtaagagtggtttgggagctatggcagcatttctgctgctgatggagcagtgaggcagtgaccaactgtgtgagagatgcaactcctgctattaagacccaaagacctaaggtccgaattcaggacgtgtaccgtgtaatgtgaaccagcttctcagctctactgtgttgacc
ctgtgtcaaggttaggggaattatagtttccaggtctccatctgtgggttaacagaggggcagagaaggcttcaaatcacatttggtaaaattgatgatgaccgcgaacaaatcgaggcccgacagtggtttgtagagagaagtcaagcaggtgactccattgtcagtcggag
tataatgcaaaagtggttgcctcaaatgtccacgggaacaaagagccttgaggtattgaaatcatgtcctagaggacatgaacgagttcaccgtgcataattctggaacggtcaatgcacacatgatgctggataaagccgtcgtgataccactgtgaaagtgtgctgtctccat
cagaagtgaatgatgctggggataacgatgagagtcaccgcaattcatccaggacgaagtgtgctgctgttgaattcatgagtgctctcatcaggaaagggaagacatcattgagctccgggtcagctctgcaccttagctgaaggccatcaagggaagcgaattgattatctgaccgttgag
ggccacatctatccactgtggtatgggagcatitaaagtctgctgattga