

HIST1H3F-histone cluster 1, H3f Gene

Reference: PTXBC069818

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H3/i; H3FI; histone H3.1; H3 histone family, member I; histone 1, H3f; histone H3/i; histone cluster 1, H3f; histone cluster 1 H3 family member f		

For research use only.

GENE DESCRIPTION	histone cluster 1, H3f
GENE ID	8968
ORIGIN SPECIES	Human
ACCESSIONS	BC069818
VECTOR	pDONR223
NCBI SYMBOL	HIST1H3F
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagcagcgaagcaaacagctcgttaagtccactggcggcaagccccgcgaagcagctggccactaaggcggctcgcaaaagcgcgcagccaccgggtggcgtgaagaagccccaccgctacaggcctggctactgtcgcctccgtgaaatccgcccgtatcagaaatcga
ctgagctactgattcgcaagctaccattccagcgtctgggtacgtgagatcgcgcaggacttcaagaccgacctgcgctccagagctcggctgtgatggcgcgtgcaggaggcctgcgaggttacctgggtgggctctttgaggacaccaacctgtgtgtatccacgccaagcgagtgac
tatcatgcccaaggacatccagctcgtcgcgcgcatcgcggagagagggcataa