



Reference: PTXBC069996

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GPR23; LPA4; P2RY9; P2Y5-LIKE; P2Y9; lysophosphatidic acid receptor 4; G-protein coupled receptor 23; LPA receptor 4; LPA-4; P2Y purinoceptor 9; P2Y5-like receptor; purinergic receptor 9		

For research use only.

GENE DESCRIPTION	lysophosphatidic acid receptor 4
GENE ID	2846
ORIGIN SPECIES	Human
ACCESSIONS	BC069996
VECTOR	pDONR223
NCBI SYMBOL	LPAR4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1