

SLC12A6-solute carrier family 12 (potassium/chloride transporters), member 6 Gene

Reference: PTXBC070107

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACCPN; KCC3A; KCC3B; solute carrier family 12 member 6; K-Cl cotransporter 3; electroneutral potassium-chloride cotransporter 3; potassium chloride cotransporter 3; potassium chloride cotransporter KCC3a-S3; potassium-chloride transporter-3a; potassium-ch		

For research use only.

GENE DESCRIPTION	solute carrier family 12 (potassium/chloride transporters), member 6
GENE ID	9990
ORIGIN SPECIES	Human
ACCESSIONS	BC070107
VECTOR	pDONR223
NCBI SYMBOL	SLC12A6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1