



Reference: PTXBC070112

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM53A; dorsal neural-tube nuclear protein; family with sequence similarity 53 member A		

For research use only.

GENE DESCRIPTION	family with sequence similarity 53, member A
GENE ID	152877
ORIGIN SPECIES	Human
ACCESSIONS	BC070112
VECTOR	pDONR223
NCBI SYMBOL	FAM53A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]