

FLJ43980-FLJ43980 protein Gene

Reference: PTXBC070117

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	FLJ43980 protein
GENE ID	124149
ORIGIN SPECIES	Human
ACCESSIONS	BC070117
VECTOR	pDONR223
NCBI SYMBOL	FLJ43980
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgaggagaagctctggctcggggccagaaggccagggccctctgctcctcaccggctgttgggtttgcagggggacgatgcgccgagagcgcgatgacagccccttcgccagtcgggtacaatatccaagagaacacctcggcaagctccagagctgccagccg
gggtgaagctccaaggtggagtgcctcttcgagcgggagcgcctgactggatgagagagacaagaagaaaaggactgctctgcatlgtgctgtgccaacgccatccagaagtggtagctctctgtggacagagagatgcaacttgatgctttgacacaaaaacaggactgct
ctgctaaggcgggtacaatgccaggagaggaatgtgcgactatttgcagaacacggagccgaccagacctccagatgtctatggcaatactacttacactatgccatcacaatgaggatcccaatgacaaaaaaactctttacacatgccaatattgaatcagcaacaagg
atgaaactcacacattttgctagctgtacatgagcaaaacagcaaatggaagacttttaagaaagcaaaaagaaaatfctcagctgaaagcattcaccaagtaatgtctgaatacaagaaaatgagacacctagaatcctcaaaatgcaatccagaaggaacatctaac
aagatggcgtcttggagaaaggagcagcaggtgcaaaagtggatgagatcctggaatcctgtaacaaggctttcaacaacacatctattgacgactcagcactatgtcagctaatgaagactttgatlttgactgaggagaagcaacagagccagcaaatgaaaaagacaaa
atggtatgggtattattgagagtgctccacaagcacaacaataatgaaacatttaa
```