

IL17F-interleukin 17F Gene

Reference: PTXBC070124

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CANDF6; IL-17F; ML-1; ML1; interleukin-17F; cytokine ML-1; interleukin 17F		
For research use only.			

GENE DESCRIPTION	interleukin 17F
GENE ID	112744
ORIGIN SPECIES	Human
ACCESSIONS	BC070124
VECTOR	pDONR223
NCBI SYMBOL	IL17F
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacagtgaagaccctgcatggcccagccatggtcaagtactgtctgtcgtatattggggcttgcctttctgagtgaggcggcagctcggaaatcccaaagtaggacatactttttccaaaagcctgagagttgccgcctgtgccaggaggtatgaagcttgacattggcatca
tcaatgaaaaccagcgcgtttccatgtcagtaacatcgagagccgtccaccccccctggaaattacactgtcactgggaccccaaccgggtacccctcggaagtgtacaggccagtgtaggaactgggctgcatcaatgctcaaggaaaggaagacatctcatgaattccgttccc
atccagcaagagaccctgggtcgtccggagggaagcaccaggctgtctgtttctttccagttggagaaggtgctgtgactgttggctgcacctgcgtcacccctgtcatccaccgtgtgcagtaa