

TMEM116-transmembrane protein 116 Gene

Reference: PTXBC070135

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 116		
For research use only.			

GENE DESCRIPTION	transmembrane protein 116
GENE ID	89894
ORIGIN SPECIES	Human
ACCESSIONS	BC070135
VECTOR	pDONR223
NCBI SYMBOL	TMEM116
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaacacaccagagtgagacagacacatctccactggatagattttactgtcagttgtcaaatggcctttgtttctcaagcctgatacctctgtattgatgacacctgtattctgtctgggaaatactagtgaatgtttccaaaactcagtcagagccacaagtgtatcttgatgcactc
accaccatcagccatggctgaactccacctctgccaacacatctgtctgtagcacattttatgtatcgcattttctgggcagctttgtactcagcctcctaccattatggcttactatccgagccagacattgtataagaattgtgaagtcactggcttctggggagtgaaca
gtgggcagtgattacattgtggaccaacgggtgcgtcttaccagtggtccttctgtctgtgggcccagctgtcattctaatgacataaagctgactaagccacaggacacaaagcttcgcatggccctttatgttctcaggctcctaacggcaacatctcagggtctactcaactgtgg
agtatatggctggacgcagcacaaattccaccaactaaagcaggagggtcggcgtgatgcagataccagacaccattattatgctcacagaagagattctatagcaggggttaaaftcactggaatccacctgacttttctgccagtactctaccatttttga