

GPR82-G protein-coupled receptor 82 Gene

Reference: PTXBC070180

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	G protein-coupled receptor 82		
For research use only.			

GENE DESCRIPTION	G protein-coupled receptor 82
GENE ID	27197
ORIGIN SPECIES	Human
ACCESSIONS	BC070180
VECTOR	pDONR223
NCBI SYMBOL	GPR82
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacaacaatacaacatgattcaaccatctatgatctctccatggcttaccatcattacatctccttctgattgtgtgttttgaacactctctctcaatggatattttaacaaaataggtaaaaaacatcaacgcacatctacctgtcacacctgtgactgcaaacctactgtgtgc
agtccatgccccttcatgagtatcttctgaaaggttccaatgggaatatcaatctgcaatgcagagtggtcaatcttctgggaactctatccatgcatgcaagatgtttgacagctcttaatttaagtggattgccataagccctatgctaccttaagcaaaagattcctcgaagag
actactcatgctatgaaaaatattatggccatttactgaaaaaattcgccagcccaactttgctagaaaactatgcaattacataggggaggtgtactggcataatcattccagttaccgtataactcagctatagaggctacagaaggagaagagagcctatgctacaatcggcagat
ggactagagccatgatctcagattgcaggtctcattggaaacacattatggatttcccttttagtagtaacatcactactcttttgaagccatctgagaaaaataagaacctgtacgtccattatggagaagattgacttacagttctgtgaaaaagacatctttgtcatccagatt
ctactaatagttgcttctctttagtattttaaacccattttatgtctacaccaaagagataactgtcagcaattgaattttaaataagaacaaaaaacattctcactgtcttctcggcagaagtagcacagacccattatattctttattagacaaaacattcaagaagacactata
atctcttcaaaagtctaattcagcacatgcaatcatatggttga