

PTP4A2-protein tyrosine phosphatase type IVA, member 2 Gene

Reference: PTXBC070182

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HH13; HH7-2; HU-PP-1; OV-1; PRL-2; PRL2; PTP4A; PTPCAAX2; ptp-IV1a; ptp-IV1b; protein tyrosine phosphatase type IVA 2; PTP(CAAXII); phosphatase of regenerating liver 2; protein tyrosine phosphatase IVA; protein tyrosine phosphatase IVA2; protein-tyrosine		

For research use only.

GENE DESCRIPTION	protein tyrosine phosphatase type IVA, member 2
GENE ID	8073
ORIGIN SPECIES	Human
ACCESSIONS	BC070182
VECTOR	pDONR223
NCBI SYMBOL	PTP4A2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaccgtccagccctgtggagatcctcatgagaacatgcgttttctgataactcacaacccctaccaatgctactctcaacaagttcacagaggaaactaagaagtatggagtacgacttgggttcgagttgtgatgctacatatgataaagctccagttgaaaaagaaggaaatccacgtt
ctagattggccatttgatgatggagctccaccctaatcagatagtagatgattgttaaacctgttaaaaaccaaatttcgtgaagagccaggttgcctgtgttcagtgcaattgtgttcaggattgggaaggccacgtgtgtgttcacttgcatttgattgaatgtggaatgaagtacgaag
atgcagttcagttataagacaaaaagaaggggagcgttcaattccaacagctgctttatttgagaaataccgacctaagatgcgattacgcttcagagataccaatgggcattgctgtgttcagtag