

IFNG-interferon, gamma Gene

Reference: PTXBC070256

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IFG; IFI; interferon gamma; immune interferon		

For research use only.

GENE DESCRIPTION	interferon, gamma
GENE ID	3458
ORIGIN SPECIES	Human
ACCESSIONS	BC070256
VECTOR	pDONR223
NCBI SYMBOL	IFNG
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaatatacaagttatcttggctttcagctctgcatcggttgggtctcttggctgttactgccaggacccatgtgaaagaagcagaaaacctaagaatattttatgcaggctcattcagatgtagcggataatggaactctttctaggcattttgaagaattgaaagaggagagtga
cagaaaaataatgcagagccaaattgtctcttttactcacaactttttaaaactttaagatgaccagagcatccaaaagagtggtggagaccatcaaggagacatgaatgtcaagttttcaatagcaacaaaaagaacgagatgacttgaaaagctgactaattatcggttaactgacttg
aatgtccaacgcaaagcaatacatgaactcatccaagtgtggtggaactgtcgcagcagctaaaacaggggaagcgaaggaagagtcagatgctgtttcgaggctgaagagcatccagtaa