



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CFHR3-complement factor H-related 3 Gene

Reference: PTXBC070259

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CFHL3; DOWN16; FHR-3; FHR3; HLF4; complement factor H-related protein 3; H factor-like 4; H factor-like protein 3; complement factor H related 3		

For research use only.

GENE DESCRIPTION	complement factor H-related 3
GENE ID	10878
ORIGIN SPECIES	Human
ACCESSIONS	BC070259
VECTOR	pDONR223
NCBI SYMBOL	CFHR3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcttactaactgcattgcacctgttgggtttctgtctataggacaagtgaaacctgtgattttccagacattaaacatggaggtctatttcatgagaatagctgtagaccatactttccagtgcctgaggaaaataattacccttattactgtgatgaacattttgagactccttcaggaa
gttactgggattacattcattgcacaaaatgggtggtcaccagcagtagaccatgtctcagaaaatgttattttcttatttggaaaattggatataacaaaatfatgggaagaagtttgcacagggttaactctacagaagttgctgccaactcctggctacggtctctcaaaagcgcagaccacagtt
acatgtacggagaagaagctggtgctctcactccagatgactcctgtgcagaacattgcctcaaaatcagatagaaaatgaaaatggattcttccgaattcctctcttatttatatttataaataagaataatataatgtaaacctggtgacacagcagatgaaatcttccacgatcaatt
catttctcagaagtgatgctcaacacaaaatttcattcactctgtattgttcattctcctcctcgcagaaaagtgttcacatctcctggggagaatgattgcttcaattctctctag

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1