

PDS5B-PDS5, regulator of cohesion maintenance, homolog B (*S. cerevisiae*) Gene

Reference: PTXBC070274

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	APRIN; AS3; CG008; sister chromatid cohesion protein PDS5 homolog B; androgen induced inhibitor of proliferation; androgen-induced proliferation inhibitor; androgen-induced prostate proliferative shutoff-associated protein AS3; androgen-induced shutoff 3;		

For research use only.

GENE DESCRIPTION	PDS5, regulator of cohesion maintenance, homolog B (<i>S. cerevisiae</i>)
GENE ID	23047
ORIGIN SPECIES	Human
ACCESSIONS	BC070274
VECTOR	pDONR223
NCBI SYMBOL	PDS5B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcattcaagactaggaccaatgatgaaaaattacatatccgctggggtcaaggaaatcagataaaatctaaaggagatggtagacgattaaagatggttgaaaaactttatggatatggaccaggactctgaagaagaaaggagctttatataacctagctttaca
tcttgcctcagatttttctcaagcatcctgataaagatgttcgcttactgtagcctgctgccttgcgtgatatttcaggatgttgcctcgaagctctacacatccctgataaaactaaaggcaagtactgatttaataactccaagattgaccgatacttggattatcttcttaa