

MKRN3-makorin ring finger protein 3 Gene

Reference: PTXBC070293

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CPPB2; D15S9; RNF63; ZFP127; ZNF127; RING finger protein 63; zinc finger protein 127; makorin ring finger protein 3		

For research use only.

GENE DESCRIPTION	makorin ring finger protein 3
GENE ID	7681
ORIGIN SPECIES	Human
ACCESSIONS	BC070293
VECTOR	pDONR223
NCBI SYMBOL	MKRN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaagagcctgcagctccctcagaagcccacgagcagccgggcccagcaggtgctgaggcagcaaggagggtgtgtctggccggacctccctgtgtgagccctccggggaatctgctgctcagattcagccctgccacatgcggcaagggctgggccccctccc
tgtagctccagtccttgcacacctccgagaggagcctgagcctgcccagcctcaggaggaggagcctggccagctccgttgccaagccgaagcagcggcatttggacaagcagatcatctgcagagatattgttga