

S100A12-S100 calcium binding protein A12 Gene

Reference: PTXBC070294

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CAAF1; CAGC; CGRP; ENRAGE; MRP-6; MRP6; protein S100-A12; EN-RAGE; calcitermin; calcium-binding protein in amniotic fluid 1; calgranulin C; extracellular newly identified RAGE-binding protein; migration inhibitory factor-related protein 6; neutrophil S100		

For research use only.

GENE DESCRIPTION	S100 calcium binding protein A12
GENE ID	6283
ORIGIN SPECIES	Human
ACCESSIONS	BC070294
VECTOR	pDONR223
NCBI SYMBOL	S100A12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacaaaacttgaagagcatctggagggaattgtcaatatctccaccaatactcagttcggaaggggcattttgacaccctctctaagggtgagctgaagcagctgcttacaaggagcttgcaaacaccatcaagaatatcaagataaagctgtcattgatgaaatattccaaggcctgg
atgctaataagatgaacaggtcgactttcaagaattcatatccctggtagccattgcgctgaaggctgccattaccacaccacaaagagtag