

XCL1-chemokine (C motif) ligand 1 Gene

Reference: PTXBC070309

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATAC; LPTN; LTN; SCM-1; SCM-1a; SCM1; SCM1A; SCYC1; SCM-1-alpha; XC chemokine ligand 1; c motif chemokine 1; chemokine (C motif) ligand 1; cytokine SCM-1; lymphotaxin; single cysteine motif 1a; small inducible cytokine subfamily C, member 1 (lymphotactin)		

For research use only.

GENE DESCRIPTION	chemokine (C motif) ligand 1
GENE ID	6375
ORIGIN SPECIES	Human
ACCESSIONS	BC070309
VECTOR	pDONR223
NCBI SYMBOL	XCL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagacttctatcctgcccctccttgccatctgctctcactgcatactgtggaaggtgtgaggagtggaagtcagataagagacctgtgtgagcctcactaccgagcgtgcccgttagcagaatcaagacctacacatcacggaaggctccttgagagcagtaattttattac
caaacgtggcctaaaagtctgtgctgafccacaagccacatgggtgagagacgtggtcaggagcatggacaggaatccaacaccagaaataacatgatccagaccaagccaacaggaaccagcaatcgaccaatcacagctgtgactctgactggttag