



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CD3D-CD3d molecule, delta (CD3-TCR complex) Gene

Reference: PTXBC070321

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD3d molecule; CD3d molecule, delta (CD3-TCR complex); CD3d antigen, delta polypeptide (TiT3 complex); CD3-DELTA; IMD19; T3D; T-cell surface glycoprotein CD3 delta chain; CD3 antigen, delta subunit; CD3 delta; OKT3, delta chain; T-cell receptor T3 delta c		

For research use only.

GENE DESCRIPTION	CD3d molecule, delta (CD3-TCR complex)
GENE ID	915
ORIGIN SPECIES	Human
ACCESSIONS	BC070321
VECTOR	pDONR223
NCBI SYMBOL	CD3D
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaacatagcacgtttctctcggcctggactggctaccctctctcgcgaagtgagccctcgaagataactatagaaagacttgaggacagagtggttggaaatgcaatacagcatcacatgggtagaggggaacgttgggaacactgctctcagacattacaagactggacctggga
aaacgacatctgcggccaccagaggaaatataatggatggagacatatatacagaagaagaatctacgcgaacatcttctcgaattgagccagagctgtgtggatcgatgacacacacacgctggcctcatctgactgattcagctgacctctgctctgtttgggagctct
ctgtctctgcacatgacatcgaagaacgtctctcggcctgcgcacacacacactctgttgggaatgacacagctctctacacgctctcagatcgcagatgctcgtcagtagccacactctgggaagaacactggcctcgaacaaactc

More info: Please contact us for any question: contact@proteogenix.fr