

AQP7P2-aquaporin 7 pseudogene 2 Gene

Reference: PTXBC070322

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	aquaporin 7 pseudogene 2
GENE ID	389756
ORIGIN SPECIES	Human
ACCESSIONS	BC070322
VECTOR	pDONR223
NCBI SYMBOL	AQP7P2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccagctgtgtctcttccatctgtggaccaggagaaacacccagcactgccaggaaacacgcactgtgataggcatcctcgtgtgtcatcatcaggggtaccatggcatgaacacaggatagccatcaatccgtcccgaggacgtcccccccgcatcttcaccttcattgtctg
gttggggcaaacgtgtcttcagggtggcatcatctacctgtgtcttcattggctccaccatccacggagagccctga