

FAM69A-family with sequence similarity 69, member A Gene

Reference: PTXBC070342

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM69A; family with sequence similarity 69 member A		
For research use only.			

GENE DESCRIPTION	family with sequence similarity 69, member A
GENE ID	388650
ORIGIN SPECIES	Human
ACCESSIONS	BC070342
VECTOR	pDONR223
NCBI SYMBOL	FAM69A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgaggagctctgtccggggcctggcctaaggaaaccctattaccctccaggctcgtctcatatgtcgggatgaaatctttctttcctggttagtggttttgggaagctggattatatgtgcagtattctacctatacagaattatgcagaggaaaggactgtaagaaataatat
gtgacaagtacaagactggagttattgatggcctgcatgtaacagcctttgtgttacagaaactcttactttggaaaatgtttatccaccaagcccaacaatcagatgtatttagggatttgggataatctaccaggtgttgaaatgtcaaatggaacaagcgttcctctgtatttggaaactga
attggaaccaagaaaagaaatgtgctatttgataagccaactagaggaaactactgtacaaaaattaaagaatggtctatagctctttaaggaccagacccggtaa