

C14orf53-chromosome 14 open reading frame 53 Gene

Reference: PTXBC070365

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf53; NCRNA00238; long intergenic non-protein coding RNA 238		
For research use only.			

GENE DESCRIPTION	chromosome 14 open reading frame 53
GENE ID	440184
ORIGIN SPECIES	Human
ACCESSIONS	BC070365
VECTOR	pDONR223
NCBI SYMBOL	C14orf53
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacaagtggcgcaactcttcaggatcttaccctcagaataagaagttctaaatagatgacaactactgaaggaaatgaatgaggacttaagctaaggaagcaggaactgctagagatgctcaaacctctagaagataaaaaaacctcttattcaaaagttaagtctaaattg
gaggaaaaacaaggagctctcagatcatgagacagatcagcagggaagggtgtgaggaaatcctcgtcatggagctccttaaggagcagagagatgaaacagaactggaaaggaaaaacaagatgctcggaaaggaaatggagatgctatggaacaagacatcgaggca
gaagaacttagtgatcaacaaaaagcaccacagacaaaaacaaggcagacttgaggatggaaggaaaaagcaacagaggaatggaatgggtcaagtatcaggacaaaaataacatcctcagaatgatttcatggcaagtgaattgagctgagaattgaagccttgaagaactac
cagaaggccaatgacctgaattatcactgtatttcagcagaattttgagccaatgcaagcatttttaaatcttctcgggtcccaagggaattagaattgtttaa