



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

HMGB3-high-mobility group box 3 Gene

Reference: PTXBC070482

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

HMG-2a; HMG-4; HMG2A; HMG4; high mobility group protein B3; high mobility group protein 2a; high-mobility group (nonhistone chromosomal) protein 4; high mobility group box 3

For research use only.

GENE DESCRIPTION	high-mobility group box 3
GENE ID	3149
ORIGIN SPECIES	Human
ACCESSIONS	BC070482
VECTOR	pDONR223
NCBI SYMBOL	HMGB3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcttaaaagggtgaccccaagaagaaaccaaaggcgcaagatgtccgtctatgctctctttgtgcagacatgcagagaagaacataagaagaaaaaccagaggctccctgtcaatttgcggaatttccaagaaggtcgtctgagaggtgggaagacgatgccggaaagagaaatctaaatttga
tgaatggcgcaaaagcagataaagtgctgatgatgcgggaatgaaggatgatggaccagctaaaggaggcagaagaagaagagatcctaattgctcccaaaaaggccacgtctggattctctgtctgtcagaattccgcccaagatcaaatccaaaaccccgcatctctatttga
gacggcgcaaaaaggctgggtgagatgtggaataaatttaattgacagtgtgaaaaacgagccttaactactaaaggctgggaagagactgaaggagaagatgaagaagagatgtgctgactataagtcgaagggaaagtgtgatgtgtgcaaaagggtcctgtcaaagttgccggaaaaagggtgga

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1