



Reference: PTXBC071583

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DAP-150; DP-150; P135; dynactin subunit 1; 150 kDa dynein-associated polypeptide; dynactin 1 (p150, glued homolog, Drosophila)		

For research use only.

GENE DESCRIPTION	dynactin 1 (p150, glued homolog, Drosophila)
GENE ID	1639
ORIGIN SPECIES	Human
ACCESSIONS	BC071583
VECTOR	pDONR223
NCBI SYMBOL	DCTN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1