

FRMPD2L2-FERM and PDZ domain containing 2 like 2 Gene

Reference: PTXBC071635

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	FERM and PDZ domain containing 2 like 2
GENE ID	728603
ORIGIN SPECIES	Human
ACCESSIONS	BC071635
VECTOR	pDONR223
NCBI SYMBOL	FRMPD2L2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacatctatcccttccaggtgaccgactcctgcaggtggatggagtgattctgtgcggcctcaccacaagcaggctgtgcagtgccctgaaggctcctgggcagggtgcaagactgggtcttagagagaagagtcgccagggtacacagcagtgctcttctgtaatgacagcatgg
gagatgaacgcacggctgttcttggtaacagccttgcctggcaggccttcgagctgtgtctcgtgacagatgggtcctaaagtgaagtaaaactaaaaagaatgccaatgggttgggattcagttctgtgagatggagaagagagctgcagccatcctcaaaagtgtctgtgagatt
aagaggtcttccgggagcagcctgagggagaatggggcattgacagctgtgacattatcctggcctgaatggaagggtccacggaaggcctcatctccaggaggtgctgacattactgagaggggccccacaggaagtcacgctcctcttggcgacccccctccaggtgcgctg
cctgagatggagcaggaatggcagacactgaactctcagctgacaaagaattcaccagggaacatgtactgactcatgtaccagccccctcctggatcaagaggacagctggaggagacagtgccctcccagatgcagggggaaggcctgggtctcaggccagagcttcccaaaag
gccatcagagagcacaatggggccaaaacagagagagaccttggccagttcctgacacattctcctgagtcaccctcattatgcaaaactaccaaagaaggatgaatcaacattggcgaccttctgaaaaaggatgtgaggcaaaactgctattcagtttgtatcatgaga
cttggaaatattctctatcctccttaaccagactttcgacagatattttctga