

SURF1-surfeit 1 Gene

Reference: PTXBC071658

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SURF1, cytochrome c oxidase assembly factor; CMT4K; surfeit locus protein 1; surfeit 1		

For research use only.

GENE DESCRIPTION	surfeit 1
GENE ID	6834
ORIGIN SPECIES	Human
ACCESSIONS	BC071658
VECTOR	pDONR223
NCBI SYMBOL	SURF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgggtggctgcttgcagctggggctcgggcggcgggctgggacgggccccggccagcggcgcctggaggagcgtcctcagggtctccccggccagggtggcctggaggccaagcagatgtggcattctgcagcagaagcatctgcccaaaagcgggaag
atgactcctttctcagtggtctcctcctcctctgactgcttggcttgggacatggcagggtccagcgtcggaaagtggaaactgaacctgattgagagctggagtcagagttctggctgagcctgtcccctgccagccgacccaatggaactgaaaaatctggagtatagc
cagtgaaagtcagggggtgcttggaccattccaaggagctgtatgatgccccggaccatgggtggaccctgctccggaggccccggggggcggcctcactcctcctcaactcagagtgggcctatgtgtcactccctccactgcacgacctgggagtcacctctgttaaatag
agggttcgtccaggaaagagtgaaactgaaacgcggcagaagggcagattgaggagaagtggaactcattgggatggtgaggctgacagaaccaggcagcctttgtccctgagaacaatccagaaggaaccactggcattatcgagactggaagctatggccagaatc
acaggcgcagagcccatctcattgatccaactccagagcacagtcctggaggaccattggaggcaaacagagttactctgaggaacgagcatctcagtagatctgacctggtatggactctcagctacatcctcctgtggttaagaaatctctactgtggacacctgt
gtgtga