

RPL15-ribosomal protein L15 Gene

Reference: PTXBC071672

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DBA12; EC45; L15; RPL10; RPLY10; RPYL10; 60S ribosomal protein L15		

For research use only.

GENE DESCRIPTION	ribosomal protein L15
GENE ID	6138
ORIGIN SPECIES	Human
ACCESSIONS	BC071672
VECTOR	pDONR223
NCBI SYMBOL	RPL15
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtgcatacaagtagatccaggagctatggagaagaagcagctctgatgtcatgcgtcttctctgaggggtccgtctgtggcagtagccagctctctctccacagggtccccgccccacccggcctgataaagcgccgactgggtctacaaggccaagcaaggttacgtt
atatataggattcgtgttcgctgtgtgcccgaacgccagttcctaaagggtgcaacttaaggcaagcctgtccatcatggtgttaaccagctaaagttgtctcgaagccttcagtcctgttcagaggagcgagctggagccactgtggggtctctgagagtcctgaattcttactgggttg
gtgaagattccacatacaaaatttttgaggttatcctcattgatccattccataaagctatcagaagaatcctgacacccagtggtatcaccacacagtcacacagcagggagatgcgtgggctgacatctgcaggccgaagagccgtggccttggaaaggccacaagtccaccac
actattgtgtgctctcgccggcagcttggagaaggcgcaataactctccagctccaccgttaccgctaa