

RPS16-ribosomal protein S16 Gene

Reference: PTXBC071674

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	S16; 40S ribosomal protein S16; ribosomal protein S16		
For research use only.			

GENE DESCRIPTION	ribosomal protein S16
GENE ID	6217
ORIGIN SPECIES	Human
ACCESSIONS	BC071674
VECTOR	pDONR223
NCBI SYMBOL	RPS16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgtccaagggctcgtcagtcgggtcaggctcttcggacgcaagaagacagcgacagctgtggcgcaactgcaaacgcggcaatggctcatcaagggtgaacggcgccctggagatgattgagccgcgcacgtacagtacaagctgtggagccagttctcttcggca
aggagcgaatttgctgtagacatccgtgtccgtgtaaggggtggtgcacgtggccagatttatgggtgagtcgccaggaactggcgcatggagaggtggctctgggagggagcccttcacagcgctctgtacccttaattgtgtgtcttttcacagctatccgtcagtcacatcc
aaagccctgggtgcctattaccagaaatgtggtgaggtccaagaaggagatcaaaagacatcctcatccagtatgaccgacccgtggtgtag