

SPIN2A-spindlin family, member 2A Gene

Reference: PTXBC071694

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SPIN2; TDRD25; dJ323P24.1; spindlin-2A; centromeric spin2 copy; spindlin homolog 2; spindlin-like protein 2; spindlin-like protein 2A; spindlin family member 2A		

For research use only.

GENE DESCRIPTION	spindlin family, member 2A
GENE ID	54466
ORIGIN SPECIES	Human
ACCESSIONS	BC071694
VECTOR	pDONR223
NCBI SYMBOL	SPIN2A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagaccccaacgcacaggaagccgaagggaacaaccagggcagctgcaggacgggccaactgggtctgcaaacatgacaaagaaaaagtctccaaagaagcagagaggcagaccttcaccaacccgcaggaaacatcgtggcctgcagaatttctcatggatgga
agggaaggagatgagcccatcacgcagtggaaggaacggtctgcatcagggtgcctataaatccctctcttctggtgaaatgatggaattgactgtgtctatggactggaacttcacagagatgaaagggtttgtctctaaattctttgcacagggtgcatcatcattatggtat
gccaaccttgcaataaccataattggcaagcagtggaacacatgtttgaggagagcagtggttctaaggatgaatggaggggagtggtcttagctcaagcacctatcatgaaagcctggtttatattacatgagaagaatcctgtctgtacatgtaccagctctagatgattataaggaa
ggtagacctccgcatatgccagaatccagtgagtcctccaacagagaggagccaggaggagttgtagatggcctaataaggtaagcatgtggaaatatacacaagaagatggctccaaaggatggcatgggtcattcacaagtggaaaccaaacctctgtgtattcatcaagttga
tgatgattccatatctatgtctacgatttgggaaaaagtcctaa