

C11orf71-chromosome 11 open reading frame 71 Gene

Reference: PTXBC071695

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C11orf71; URLC7; up-regulated in lung cancer 7; chromosome 11 open reading frame 71		

For research use only.

GENE DESCRIPTION	chromosome 11 open reading frame 71
GENE ID	54494
ORIGIN SPECIES	Human
ACCESSIONS	BC071695
VECTOR	pDONR223
NCBI SYMBOL	C11orf71
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccctgaacaatgtgtccctgtcctcgggtgatcagaggagcaggggtggcctaccgtcttccatggcgacctcagaccgcccgtcagcgttggcgatgtgtcctcgggagacgcttctcgtttcaggcctgaagcgattcatctaggacctcggcaggcgggtgcgaccaagc
gttcgggcccagagccgtcagtgatgtgtggcggccggagcccaagggaaccagatggccggggcggagccgccaagccagattctcacttaccctccctgcccgttgaaccggtatctcctaagaagtgtgtcgaacagcgtttgattgcattaggagttactacgtggacct
gaagctatctgagctcgaacccctgatttaagatgggtgtcaacatgcagccattcaaatattttgaaatttttagtaatatga