



Reference: PTXBC071729

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BTBD15; HSPC063; ZNF851; zinc finger and BTB domain-containing protein 44; BTB (POZ) domain containing 15; BTB/POZ domain-containing protein 15; zinc finger protein 851; zinc finger and BTB domain containing 44		

For research use only.

GENE DESCRIPTION	zinc finger and BTB domain containing 44
GENE ID	29068
ORIGIN SPECIES	Human
ACCESSIONS	BC071729
VECTOR	pDONR223
NCBI SYMBOL	ZBTB44
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]