

FAM71E1-family with sequence similarity 71, member E1 Gene

Reference: PTXBC071737

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM71E1; family with sequence similarity 71 member E1		

For research use only.

GENE DESCRIPTION	family with sequence similarity 71, member E1
GENE ID	112703
ORIGIN SPECIES	Human
ACCESSIONS	BC071737
VECTOR	pDONR223
NCBI SYMBOL	FAM71E1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggccccccttggcctgatctccaggagcccccctccaggagctcatcacaatcaggagccgctcttggctgacgtcataaagcccctccgcatcatgacgtcacagtgcgcgtagtcgcccccctcgtttctccctctgctcctccgtccgctccgacgggg
acattgcaatgagcgggatcgccctaaagccgcccctgggtggagctggcggagtggaaccagggtggatggcagcctctcccaagggccgtccagacggctccaacgtacctccagagcggcgaatcgaccagtttcgggacttcccctcttggagcaactctgata
ggtagctcgggtgggagaagtgccaacgaggtcaccatgggggtggcagcctccagtcagccctggagctccggacattgctctggccggccctgccaaggagaacggacactgcaactctcgggtgttcccccttgaagttcgtccaagctcttgcacgacaaaagccg
gtgtcagctcgaggtcaagttgaaccaccagccgacctctacttgcagctcggccccactcaagaccgagaccgagagttggccagttgggtgctgctctaccgctcgttctctctgcttctgctgtgcccttcacgcaggagtaa