

TRGV3-T cell receptor gamma variable 3 Gene

Reference: PTXBC071739

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	T cell receptor gamma variable 3
GENE ID	6976
ORIGIN SPECIES	Human
ACCESSIONS	BC071739
VECTOR	pDONR223
NCBI SYMBOL	TRGV3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgggtggccctactgggtgctttagctttctgtctcctgcagtcagaaatctccaacttgaaggagaacgaagtcagtcaccaggcagactgggtcatctgctgaaatcaacttgcgatcttactgtaacaaatccttctacatccactggtacacaccaggagggaagggccc
cacagcgtcttctgtactatgacgtctccactgcaaggatggttgaatcaggactcagtcaggaaagtattatactacacaccaggagggtggagctggatattgagactgcaaaatctaattgaaatgattctgggtctattactgtgccacctgggacaggcagatgattcagac
ctgtctacaccacactgaaatctgcttggctgctctgtgacacaagatagagccgccctctcatttctgccaccaatttctgtattctgaacaagagaaagacagcttaactcctgatctccctcctaatacacactgtcctggcagcagctcatcctgttccccaccctccc
ccacaacttctgaagatcaagctgcatctccagcctcagctaa