

CD40LG-CD40 ligand Gene

Reference: PTXBC071754

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD154; CD40L; HIGM1; IGM; IMD3; T-BAM; TNFSF5; TRAP; gp39; hCD40L; CD40 antigen ligand; CD40-L; T-B cell-activating molecule; T-cell antigen Gp39; TNF-related activation protein; tumor necrosis factor (ligand) superfamily member 5		

For research use only.

GENE DESCRIPTION	CD40 ligand
GENE ID	959
ORIGIN SPECIES	Human
ACCESSIONS	BC071754
VECTOR	pDONR223
NCBI SYMBOL	CD40LG
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatcgaacatacaacaaactctccccgatctcgccgactggactgccatcagcatgaaaattttatgtattactactgtttttctatcaccagatgattgggtcagcacttttctgtgtatctcatagaaggttgacaaagatagaagatgaaaggaattctcatgaagatttga
ttcatgaaaacgatacagagatgcaacacaggagaaagatcttcttactgaactgtgaggagattaaaagccagtttgaaggctttgtgaaggatataatgttaacaaagaggagacgaagaagaaaacagctttgaaatgcaaaaagggtgatcagaatcctcaattgcggcacat
gtcataaagtgaaggccagcagtaaaacaacatctgtgttacagtggtgtaaaagatactacacatgagcaacaactgttaacctggaaaatgggaacagctgacctttaaagacaaggactctattatctatgcccagtcacctctgttccaatcgggaagcttcgagtc
gctccatttatagccagcctctgcctaaagtcctccggttagattcgagagaatcttactcagagctgcaaatccacagttccgcaaaccttgcgggcaacaatccattcacttggaggagatttgaattgcaaccaggtgcttcggtgttgcattgtgactgatccaagccaagtga
ccatggcactggcttcacgtcttggcttactcaaacctga