

### LYRM5-LYR motif containing 5 Gene

Reference: PTXBC071769

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LYRM5; electron transfer flavoprotein regulatory factor 1; LYR motif containing 5; LYR motif-containing protein 5		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	LYR motif containing 5
<b>GENE ID</b>	144363
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC071769
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	LYRM5
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

**ORF (Only for reference, product provided without STOP codon )**

atggccaattcttaagaggagaagtactaaaactttataaaaactctgctgtatcttgacgagactatccaaaaggagcagactatttaaaaagcgtttgaagaacatttcccttaaaaacaagatgtgaagaatccagagaagatcaagaactattgcacagggcgaatttgaatgaaagagctagaagctttgacttcccttaggaaatcacagagctatgaaacaacgctattattcagataccaaaaaactaattga