

TMEM207-transmembrane protein 207 Gene

Reference: PTXBC071780

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UNQ846; transmembrane protein 207; SRSR846		
For research use only.			

GENE DESCRIPTION	transmembrane protein 207
GENE ID	131920
ORIGIN SPECIES	Human
ACCESSIONS	BC071780
VECTOR	pDONR223
NCBI SYMBOL	TMEM207
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcaagatccagacttttcagtgctcactcagcgaatcgaacgatagggatctgtgttgcgcgtattccagttgggtgctcggacctaccatgcgaagaagatgaaatgtgtgtaattataatgaccaacacccctaatggctgtatctggtatcctctgctgctggtttgtggcagct
cttctctgtggagctgtgtcctctgctcctcagtgctggctgaggagacccgaattgattctcacaggcgcaccatggcagtttttctgttggagacttggactctattatgggacagaaagcagctgtgagtcacactgttggaaatcaccttcaaaactcaaacccctgacctatctgttcc
tgctcatgttttggcccttaggctccccacccctccatgaagaaattgtaaaaacaacctga