

NCOA7-nuclear receptor coactivator 7 Gene

Reference: PTXBC071782

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NCOA7-AS; ERAP140; ESNA1; Nbla00052; Nbla10993; TLDC4; dJ187J11.3; nuclear receptor coactivator 7; 140 kDa estrogen receptor-associated protein; TBC/LysM-associated domain containing 4; estrogen nuclear receptor coactivator 1; estrogen receptor associated		

For research use only.

GENE DESCRIPTION	nuclear receptor coactivator 7
GENE ID	135112
ORIGIN SPECIES	Human
ACCESSIONS	BC071782
VECTOR	pDONR223
NCBI SYMBOL	NCOA7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggataccaaggaagagaaggaacggaacaaagtattttgctcgactgaaaaagaaaaacaagccaaacaaatgcagagacagcctcagctgtagctacaaggactcactatgggaaggaagataataacagtagtttagagccagacaagtgcacattgctgtgga
agaggaatataatgactgatgagaaaaaaagagaaaaagtaatcagttaaaggagatcaggcgtacagaactaaagagatattatagtagttggtgagttatcatggcgttactgcagtattttcatttaa