

TNS3-tensin 3 Gene

Reference: PTXBC071791

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TEM6; TENS1; tensin-3; tensin-like SH2 domain containing 1; tensin-like SH2 domain-containing protein 1; thyroid specific PTB domain protein; tumor endothelial marker 6; tensin 3		

For research use only.

GENE DESCRIPTION	tensin 3
GENE ID	64759
ORIGIN SPECIES	Human
ACCESSIONS	BC071791
VECTOR	pDONR223
NCBI SYMBOL	TNS3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggaggccatgggtggacacactacatcacgagcgcacatcgtgtgtccctccctgccgctgctctgaggagtcctacctgcacacctacagagggtcacgcgcgtcgaagccaagcaggggacaactacctggtattaaaccttcagaaaagagatatgacc
ttacgaagcttaacccaagatcatggatgtgggtggccagagctccacgcaccgcccctggataagatgtgtacatatacgaagcgaggagtcctggctgaacagcaacctccagcatgtgtgtcgtcattcactcagggggcgggaaaggacgcataggagtggtcatatcatct
acatgcatttcaccaacgtctcagccagcggcaccagcccttgacaggtttgcaatgaagaagtttatgatgacaaagtttcagctttaaagcagcctcccaaaacggatgttcagttcctcagtggtgctcctgtccggatcgggtgaaaatgaatgcctctccctgttctcgtatgttc
atctccacggcaccgccaaactcgcacaggtggaggtgtgcccggcccttttgaagctctaccaagccatgcagcctgtgtacacctccggatcacaacgttgcccagaaaacccagcaggtatcgtatcgtcagccggccagctctgaaggagatgtcatgtcatga
attacaatttgcaacatctga