

CCDC75-coiled-coil domain containing 75 Gene

Reference: PTXBC071798

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CCDC75; CENP-Y; CENPY; G patch domain-containing protein 11; centromere protein Y; coiled-coil domain containing 75; coiled-coil domain-containing protein 75; G-patch domain containing 11		

For research use only.

GENE DESCRIPTION	coiled-coil domain containing 75
GENE ID	253635
ORIGIN SPECIES	Human
ACCESSIONS	BC071798
VECTOR	pDONR223
NCBI SYMBOL	CCDC75
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccaggaagtgccattggtcatgaggcatcattaaacggaaagcagaggaaaaattgaaagctacagaaaaagattcacatgaaaaaccaagctgaagaaaaagctgcagaacgttcgaatgcgactaaaaataagcaagatgaaatgaagctagaaggagatctcag
aagaagccagcgagcctgtcaacaactggatgtccagaaaaatattcaggttccaggggaagcatggtactggttgaagcctgaagaggagactgaagaagatgaagaagaaaaagacaggatgaagatgaatataagagtggaagatttaagcgtactggaaaaattacaataattgac
tagttatttaagagaagaacatctatattgtattgtgtggaacgcctatgaagataaagaacacatcttcaaatgtccaggaccaacttctgcagatcatgactaa