

HOXA10-homeobox A10 Gene

Reference: PTXBC071843

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	homeobox protein HOXA10; HOX1; HOX1.8; HOX1H; homeobox protein Hox-A10; homeo box A10; homeobox protein 1H; homeobox protein Hox-1.8; homeobox protein Hox-1H; homeobox A10		

For research use only.

GENE DESCRIPTION	homeobox A10
GENE ID	3206
ORIGIN SPECIES	Human
ACCESSIONS	BC071843
VECTOR	pDONR223
NCBI SYMBOL	HOXA10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtcaaggcaattccaagggtgaaacgcagccaactggctcagggcaagagtggtcgggaagaagcgtgcccctacacgaagcaccagactggagctggagaaggagttctgttcaatgtacctactcgagagcggcgctagagattagccgcagcgtccacctcac
ggacagacaagtgaatactggttcagaaccgcagatgaactgaagaaatgaatcgagaaaaccgatccggagctcacgccaacttaattttctga