



Reference: PTXBC071859

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	large neutral amino acids transporter small subunit 4; L-type amino acid transporter 4; solute carrier family 43 (amino acid system L transporter), member 2; solute carrier family 43 member 2		

For research use only.

GENE DESCRIPTION	solute carrier family 43, member 2
GENE ID	124935
ORIGIN SPECIES	Human
ACCESSIONS	BC071859
VECTOR	pDONR223
NCBI SYMBOL	SLC43A2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]