



Reference: PTXBC071865

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LOXL; lysyl oxidase homolog 3; lysyl oxidase-like protein 3; lysyl oxidase like 3		
For research use only.			

GENE DESCRIPTION	lysyl oxidase-like 3
GENE ID	84695
ORIGIN SPECIES	Human
ACCESSIONS	BC071865
VECTOR	pDONR223
NCBI SYMBOL	LOXL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr