



Reference: PTXBC071906

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CK5; DDD; DDD1; EBS2; KRT5A; keratin, type II cytoskeletal 5; 58 kda cytokeratin; CK-5; cytokeratin-5; epidermolysis bullosa simplex 2 Dowling-Meara/Kobner/Weber-Cockayne types; keratin 5 (epidermolysis bullosa simplex, Dowling-Meara/Kobner/Weber-Cockayne		

For research use only.

GENE DESCRIPTION	keratin 5
GENE ID	3852
ORIGIN SPECIES	Human
ACCESSIONS	BC071906
VECTOR	pDONR223
NCBI SYMBOL	KRT5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr