

MARK4-MAP/microtubule affinity-regulating kinase 4 Gene

Reference: PTXBC071948

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MARK4 serine/threonine protein kinase; MARK4L; MARK4S; MARKL1; MARKL1L; PAR-1D; MAP/microtubule affinity-regulating kinase 4; MAP/microtubule affinity-regulating kinase like 1; microtubule affinity regulating kinase 4		

For research use only.

GENE DESCRIPTION	MAP/microtubule affinity-regulating kinase 4
GENE ID	57787
ORIGIN SPECIES	Human
ACCESSIONS	BC071948
VECTOR	pDONR223
NCBI SYMBOL	MARK4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcttcgacgggtgctgccccgggcaacgacgcgaactcggacacggctgagaacctctgctggatgccgagggccaacatcaagattgctgacttggcttcagcaacgagttcacgctgggatcgaagctggacagttctgcccggagcccccatatgccgccccggagctgtttcagggaagaagtacgacggccggagggtggacatctggagcctgggagtcacctgtacacctcgtcagcgctccctgcccttcacgggcacaaacctcaagaggtga