



Reference: PTXBC071958

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NET49; NOC4; UTP19; nucleolar complex protein 4 homolog; NOC4 protein homolog; NOC4-like protein; nucleolar complex-associated protein 4-like protein; nucleolar complex associated 4 homolog		

For research use only.

GENE DESCRIPTION	nucleolar complex associated 4 homolog (S. cerevisiae)
GENE ID	79050
ORIGIN SPECIES	Human
ACCESSIONS	BC071958
VECTOR	pDONR223
NCBI SYMBOL	NOC4L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]