

OAS1-2',5'-oligoadenylate synthetase 1, 40/46kDa Gene

Reference: PTXBC071981

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E18/E16; IFI-4; OIAS; OIASI; 2'-5'-oligoadenylate synthase 1; (2-5')oligo(A) synthase 1; 2',5'-oligo A synthetase 1; 2'-5'-oligoadenylate synthetase 1, 40/46kDa; 2'-5'-oligoisoadenylate synthetase 1; 2-5A synthase 1; 2-5A synthetase 1; E16 (2'-5') oligo A		

For research use only.

GENE DESCRIPTION	2',5'-oligoadenylate synthetase 1, 40/46kDa
GENE ID	4938
ORIGIN SPECIES	Human
ACCESSIONS	BC071981
VECTOR	pDONR223
NCBI SYMBOL	OAS1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatggatctcagaaataccacagccaaatctctggacaagttcattgaagactatctctgccagacagctgtttccgcatgcaaatcaaccatgccattgacatcatctgtgggttcctgaaggaaaggcttcagggtagctcctaccctgtgtgtgtccaaggtgtaagggtgg
ctctcaggcaaggccaccacctcagaggccgactcagcgtgacctggtgtctctcagtcctctcaccactttcaggatcagttaaatcgccgggagagttcatccaggaaattaggagacagctggaagcctgtcaaaagcagagagcatttccgtgaagtgtgaggtccagg
ctccagctggggcaacccccgtgcctcagctctgactgagttcgtccagctcggggaggggggtggagttcgtatgtctgctccttctgagccctgggtcagttgactggcagctataaacctaaccaccaatctatgcaagctatcaggagtgaccgacctgcagaaaga
ggcgagttctccacctgcttcacagaactacagagagactctgaagcagcggccccaagctcaagagcctatccgctagtcgaagcactggtacccaaattgaaagaagcgtgggaagctgccacctcagatgccctggagctctgacggtctatggtggagcgagg
gagcatgaaaacacatttcaacacagccagggaattcggagcgtcttggaaattagtcataaactaccagcaactctgcatctactggacaagattatgactttaaaaccccatattgaaaagtacctgagaaggcagctcacgaaacccaggcctgtgatctggacccgggaccc
tacaggaaactgggtggaggagacccaagggttggaggcagctggcacaagaggctgagggcctggctgaattaccatgctttaagaattggagtggtgccagtgagctcctggaattctgctggtgagacctcctgtctcctcctgccattcatccctgccctctcatgaagcttg
a